



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hybrid Electro-Hydraulic and Electro-Pneumatic Equipment

**Product Code :**  
LBNY-0006-2110002



**Description :**

Hybrid Electro-Hydraulic and Electro-Pneumatic Equipment Teaching Equipment Mechatronics Trainin

**Technical Specification :**

The experimental device uses double-sided structure, ie an experimental table provides places for doing pneumatic and hydraulic test experiment for two groups of four students. It can simultaneously do double-sided comprehensive pneumatic and hydraulic experiments, optimizing and sharing resources and improving products' performance. Its overall structure is simple and practical.

The training equipment is iron double matte dense spray pattern, which is beautiful, durable and practical. The basic configuration is as following: 1 set of pneumatic training components; 1 set of pneumatic PLC control module; 1 piece of low noise air compressor; 1 set of hydraulic training components; 1 set of hydraulic PLC control module; 1 piece of hydraulic test pump. The optional configuration is the computer.



**Laboratory Instrument India**