



Email : sales@laboratoryinstrumentindia.com

Product Name :
Proportional Pneumatic Trainer

Product Code :
LBNY-0006-2120001



Description :

Proportional Pneumatic Trainer, Proportional Pneumatic Trainer Equipments, Proportional Pneumatic Trainer Tools, Proportional Pneumatic Trainer Tool Kits, Proportional Pneumatic Trainer Manufacturers, Proportional Pneumatic Trainer Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Gabon, Lesotho, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia

Technical Specification :

The design of double-sided experiment platform can satisfy two groups of students at the same time. The components, circuits and electrical control of the two sides are independently controlled, and the experiment is not interfered with each other.

The modular structure design of structures, the experiment is simple, convenient, various pneumatic components into separate modules, with convenient installation plate, the experiment can freely on general aluminum plate, when assembling various experimental circuits, the operation is simple and convenient. Quick and reliable connection, easy installation and convenient connection.

Standard industrial components, reliable edge and safety.

The laboratory has a computer communication interface which can be controlled with PLC programming device; The gas source adopts the no-oil silent air compressor, which has the characteristics of low noise (Multiple

control methods: single pneumatic control circuit, electric control pneumatic circuit with relay, PLC control pneumatic circuit, etc.



Laboratory Instrument India