



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Pneumatic and Electro-Pneumatic Bench

**Product Code :**  
LBNY-0006-2120002



**Description :**

Pneumatic and Electro-Pneumatic Bench School Teaching

**Technical Specification :**

Pneumatic and electro-pneumatic bench Double-sided mobile workstation allows the assembly and operation of the various components without the use of tools (components equipped with quick assembly system)

1 Rolling table dimension between length 1200 mm, height 1700 mm and depth 800 mm minimum.

1 extruded and anodized aluminum plate of dimension between 1080 mm length and 680 mm depth minimum with the possibility of fixing the electrical components on both sides.

2 cabinet with a minimum of four drawers and locks

1 cable guide

1 x power supply unit for mounting frame 220V AC- 24 VDC/ 4A minimum with minimum 6 outputs

1 x silent lubricated compressor for classroom use: protected by pressure switch and safety valve



**Laboratory Instrument India**