



Email : sales@laboratoryinstrumentindia.com

Product Name :
PIC Trainer

Product Code :
LBNY-0006-2140002



Description :

PIC Trainer, Didactic Equipment, Microprocessor Trainer Kit

Technical Specification :

It's mainly suitable for learning of PIC SCM. equipped with three kinds of PIC16F877A, PIC16F876A, PIC16F628A chip minimum system, normal input output experiment unit, kit self equipped programmer module, convenient for chip burning program. SCM minimum system burning interface supports connect simulator to do simulation.



Laboratory Instrument India