



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Oil Volumetric Pump Study Bench

**Product Code :**  
LBNY-0006-2150002



**Description :**

Oil Volumetric Pump Study Bench Teaching Aid Equipment Heat

**Technical Specification :**

Aluminum profile frame structure, bottom with universal wheel, flexible movement. There are three pumps, one piston pump, one gear pump, one rotating blade pump, oil storage tank and piping system.

Experimental contents:

The characteristic test and operation of the three pumps are compared, the circuit pressure, temperature and flow measurement of different pumps are determined, and the electrical parameters of different pumps are determined.

The bench is equipped by 3 pumps: one piston pump, one gear pump and a rotary vane pump. It allows the study of the characteristics and the performances of the volumetric pumps.



**Laboratory Instrument India**