



Email : sales@laboratoryinstrumentindia.com

Product Name :
Mechanical Alignment

Product Code :
LBNY-0006-2160001



Description :

Mechanical Alignment Technology Integrated Trainer Didactic, Mechanical Alignment Equipments, Mechanical Alignment Tools, Mechanical Alignment Tool Kits, Mechanical Alignment Manufacturers, Mechanical Alignment Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Gabon, Lesotho, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia

Technical Specification :

The device is in accordance with relevant national occupational standards and industry standards, combined with various vocational schools, technical colleges, "CNC technology and its application", "machinery manufacturing technology", "mechanical and electrical equipment installation and maintenance," "mechanical assembly", "mechanical equipment assembly and automatic control "and other professional training objectives developed. Mainly trains the student to read and draw the assembly drawing and part drawing, fitter basic operation, spare parts and the assemble process and adjustment, assembling quality inspection skills. Improve the employability of the students in the machinery manufacturing enterprises and related industries assembly line process and implementation, installation, commissioning and maintenance of electrical and mechanical

equipment , machining quality analysis and control, primary production and management positions.



Laboratory Instrument India