



Email : sales@laboratoryinstrumentindia.com

Product Name :
Digital IP CCTV

Product Code :
LBNY-0002-670005



Description :

Digital IP CCTV

Technical Specification :

Designed to withstand all elements and is highly damage-resistant.
Available with an additional privacy feature that adds to preventing any backdoor entries.
Can be used both indoors and outdoors without any trouble anywhere, Comes in a discrete design that makes it difficult to determine where the camera is pointing.
1080 P IP CAMERA

Package includes:

1x4Channel H-264 Network Video Recorder(POENVR)
1x1080(2 Megapixels) IR Dome Day /Night IP CAMERA
1x1000 GB SATA Hard Disk
Cable and Accessories.

Digital IP CCTV, Digital IP CCTV Equipments, Digital IP CCTV Tools, Digital IP CCTV Tool Kits, Digital IP CCTV Manufacturers, Digital IP CCTV Suppliers from India, China, Kenya

Laboratory Instrument India

