



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sanitation Fittings Training Panel

Product Code :
LBNY-0006-2170001



Description :

Sanitation Fittings Training Panel Teaching Aid Equipment Fluids, Sanitation Fittings Training Panel Equipments, Sanitation Fittings Training Panel Tools, Sanitation Fittings Training Panel Tool Kits, Sanitation Fittings Training Panel Manufacturers, Sanitation Fittings Training Panel Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Gabon, Lesotho, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia

Technical Specification :

Sanitation fittings can be connected and investigated in respect of their function and operation on this trainer. The trainer provides several methods of connection so additional products from different manufacturers can be used. Some connectors are flexible and will adapt to different sized fittings.

Technical details

Specification

investigating the most common sanitation fittings
facilities for hot and cold water supply using quick-release couplings
possibility to test different fittings simultaneously
flushing valve

two handle mixer for a shower
two handle mixer as water tap
3 dial thermometers
2 Bourdon tube manometers
2 rotameters
2 adjustable pressure reducing valves
Technical data
Flushing valve



Laboratory Instrument India