



Email : sales@laboratoryinstrumentindia.com

Product Name :
Model Sedimentation Tank

Product Code :
LBNY-0006-2180002



Description :

Model Sedimentation Tank Vocational Training Equipment

Technical Specification :

This equipment has been designed to study the hydraulic behaviour and the settling efficiency of a sedimentation tank.
Framework of AISI 304 stainless steel with castors
Tank of clear methacrylate with capacity of 80 litres and section for mixing tap water and suspension
Suspension feed tank with capacity of 120 litres and recycle
Centrifugal pump of AISI 304 stainless steel
Variable area flowmeter for tap water
Variable area flowmeter for suspension
Power supply: 230 Vca 50 Hz single-phase - 0.5 KVA
(other voltage and frequency on demand)



Laboratory Instrument India

