



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Model Sedimentation Tank

**Product Code :**  
LBNY-0006-2180002



**Description :**

Model Sedimentation Tank Vocational Training Equipment

**Technical Specification :**

This equipment has been designed to study the hydraulic behaviour and the settling efficiency of a sedimentation tank.  
Framework of AISI 304 stainless steel with castors  
Tank of clear methacrylate with capacity of 80 litres and section for mixing tap water and suspension  
Suspension feed tank with capacity of 120 litres and recycle  
Centrifugal pump of AISI 304 stainless steel  
Variable area flowmeter for tap water  
Variable area flowmeter for suspension  
Power supply: 230 Vca 50 Hz single-phase - 0.5 KVA  
(other voltage and frequency on demand)



**Laboratory Instrument India**

