



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Sorting and Packaging Machine Trainer

**Product Code :**  
LBNY-0006-2200002



**Description :**

Sorting and Packaging Machine Trainer, Teaching Equipment

**Technical Specification :**

**Brief Introduction**

automatic granule packing machine is mainly used packing different kinds of coarse powder?granule material side sealing flat bag or pillow bag filling and packing. Specific materials that can be packed: all kinds of soup powder, seasonng, Chinese patent medicine, coffee, soybean, solt, melon seeds, tea, etc..

Products Composition and Packing Process

General Composition



**Laboratory Instrument India**