



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sorting and Packaging Machine Trainer

Product Code :
LBNY-0006-2200002



Description :

Sorting and Packaging Machine Trainer, Teaching Equipment

Technical Specification :

Brief Introduction

automatic granule packing machine is mainly used packing different kinds of coarse powder?granule material side sealing flat bag or pillow bag filling and packing. Specific materials that can be packed: all kinds of soup powder, seasonng, Chinese patent medicine, coffee, soybean, solt, melon seeds, tea, etc..

Products Composition and Packing Process

General Composition



Laboratory Instrument India