



Email : sales@laboratoryinstrumentindia.com

Product Name :
Inverter (3000 W)

Product Code :
LBNY-0002-670006



Description :

Inverter (3000 W)

Technical Specification :

With charging function that could charge your battery
Makes a best solution for urgent power supply demand
With Overload protection and overheat protection
Voltage: 220V, 50 Hz. Output
Input 12 VDC.

Inverter (3000 W), Inverter (3000 W) Equipments, Inverter (3000 W) Tools, Inverter (3000 W) Tool Kits, Inverter (3000 W) Manufacturers, Inverter (3000 W) Suppliers from India, China, Kenya



Laboratory Instrument India