



Email : sales@laboratoryinstrumentindia.com

Product Name :
UHT HTST System

Product Code :
LBNY-0005-1360001



Description :

UHT HTST System, UHT HTST System Equipments, UHT HTST System Tools, UHT HTST System Tool Kits, UHT HTST System Manufacturers, UHT HTST System Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Gabon, Lesotho, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia

Technical Specification :

Touch screen full colour control panel (High resolution widescreen (800 x 480px))
Graphical display of temperatures
32 storage memories (presets)
Temperature data can be saved to USB data stick (as standard)
Advanced data logging and display option for saving to PC
Calculated flow display (as standard)
Electromagnetic flow sensor option
Hold time display
F0 display



Laboratory Instrument India